

Product datasheet

Anti-17 beta Estradiol antibody ab50660

Overview

Product name	Anti-17 beta Estradiol antibody
Description	Goat polyclonal to 17 beta Estradiol
Host species	Goat
Specificity	Cross reactivities with related molecules; Estrone < 1%, Estriol < 1%, Testosterone < 1%, DHEA < 1%, Androstenedione 0%.
Tested applications	Suitable for: RIA, ELISA
Immunogen	17 β Estradiol-3 conjugated to BSA

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
Storage buffer	Preservative: 0.1% Sodium azide Constituent: Whole serum
Purity	Whole antiserum
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab50660** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
RIA		
ELISA		

Application notes	EIA: Titre 1/1100000. RIA: Titre 1/1100000.
--------------------------	--

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

Target

Relevance	17 beta Estradiol is a steroid hormone that is produced in the ovary, placenta, testis, and possibly the adrenal cortex.
Cellular localization	Secreted

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors