

Product datasheet

Anti-Biotin antibody (FITC) ab53469

★★★★☆ 1 Abreviews 1 Image

Overview

Product name	Anti-Biotin antibody (FITC)
Description	Rabbit polyclonal to Biotin (FITC)
Host species	Rabbit
Conjugation	FITC. Ex: 493nm, Em: 528nm
Specificity	Reacts with free biotin and biotinylated proteins regardless of species.
Tested applications	Suitable for: ICC/IF
Immunogen	d-Biotin conjugated to KLH
General notes	Fluorescein isothionate; Ab/FITC ratio=approximately 1/20.

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	Preservative: None Constituents: 50% Glycerol, PBS
Purity	Protein A purified
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab53469** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


Application	Abreviews	Notes
ICC/IF	★★★★☆	1/100.

Target

Relevance

Biotin is a water soluble vitamin, generally classified as a B complex vitamin, also called vitamin B4. After the initial discovery of biotin, nearly forty years of research were required to establish it as a vitamin. Biotin is required by all organisms but can only be synthesized by bacteria, yeasts, molds, algae, and some plant species. Biotin is required as prosthetic group of enzymes involved in incorporation of carbon dioxide into organic compounds. Biotin has a MW of 244 Da.

Images



ab53469 staining Biotin in human ovarian cells by Immunocytochemistry/ Immunofluorescence. Cells were fixed with formaldehyde, permeabilized with 0.1% Triton $\times 100$ in PBS and blocking with 5% serum for 30min at 23°C was performed. Samples were incubated with primary antibody (1/100) for 1 hour at 23°C .

Immunocytochemistry/ Immunofluorescence - Anti-Biotin antibody (FITC) (ab53469)

This image is a courtesy of Anonymous Abreview

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors