

## Product datasheet

# Anti-C7ORF61 antibody ab84941

1 Image

### Overview

<b>Product name</b>	Anti-C7ORF61 antibody
<b>Description</b>	Rabbit polyclonal to C7ORF61
<b>Host species</b>	Rabbit
<b>Tested applications</b>	<b>Suitable for:</b> WB
<b>Species reactivity</b>	<b>Reacts with:</b> Human
<b>Immunogen</b>	Synthetic peptide within residues: YLVWAVRKH LRRLYRRQER HRRHHV RCHA APRPNPAQSL KLDAQSPL, corresponding to C terminal amino acids 159-208 of Human C7ORF61 (NP_001004323) <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>
<b>Positive control</b>	Human fetal brain lysate.

### Properties

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
<b>Storage buffer</b>	Preservative: 0.09% Sodium azide Constituents: 2% Sucrose, PBS
<b>Purity</b>	Immunogen affinity purified
<b>Purification notes</b>	ab84941 is purified by a peptide affinity chromatography method.
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG

### Applications

Our [Abpromise guarantee](#) covers the use of **ab84941** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		Use a concentration of 1 µg/ml. Detects a band of approximately 31 kDa (predicted molecular weight: 24 kDa). Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

## Target


---

### Relevance

The exact function of C7ORF61 remains unknown.

## Images

---



Anti-C7ORF61 antibody (ab84941) at 1  $\mu$ g/ml  
(in 5% skim milk / PBS buffer) + Human fetal  
brain lysate at 10  $\mu$ g

### Secondary

anti-Rabbit IgG HRP at 1/50000 dilution

**Predicted band size:** 24 kDa

**Observed band size:** 31 kDa

Gel concentration 12%

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors