

Product datasheet

Anti-CCNO antibody ab47681

★★★★☆ 1 Abreviews 1 Image

Overview

Product name	Anti-CCNO antibody
Description	Rabbit polyclonal to CCNO
Host species	Rabbit
Specificity	At least two isoforms of Cyclin O are known to exist; this antibody will recognize both isoforms.
Tested applications	Suitable for: WB
Species reactivity	Reacts with: Mouse, Rat, Human
Immunogen	A 13 amino acid peptide from near the center of Human Cyclin O (Genbank accession No. NP_066970).
Positive control	Human bladder tissue lysate.

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	Preservative: 0.02% Sodium azide Constituent: PBS
Purity	Immunogen affinity purified
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab47681** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		Use a concentration of 1 - 2 µg/ml. Detects a band of approximately 40 kDa (predicted molecular weight: 38 kDa).

Target


Sequence similarities

Belongs to the cyclin family.

Developmental stage

Maximum levels during G(1) phase. Levels decrease through S and G(2) phases.

Images



Western blot analysis of CCNO in Human bladder tissue lysate with [ab134082](#) at (1) 1 and (2) 2 $\mu\text{g/mL}$.

Western blot - Anti-CCNO antibody (ab47681)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors