

Product datasheet

Anti-COX1 / Cyclooxygenase 1 antibody ab39351

Overview

Product name	Anti-COX1 / Cyclooxygenase 1 antibody
Description	Rabbit polyclonal to COX1 / Cyclooxygenase 1
Host species	Rabbit
Specificity	This antibody is specific for Prostaglandin I synthase.
Tested applications	Suitable for: WB, IP
Species reactivity	Reacts with: Sheep, Cow, Human Does not react with: Rat
Immunogen	Synthetic peptide: LLKNPEALAA VRGELETVLL GAEQPISQMT T conjugated to KLH, corresponding to amino acids 299/329 of Cow Prostaglandin I synthase Run BLAST with ExPASy Run BLAST with NCBI

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Storage buffer	Preservative: 0.02% Sodium Azide Constituents: 50% Glycerol, Tris buffered saline, 0.5mg/ml BSA, pH 8.0
Purification notes	ab39351 is an affinity purified IgG.
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab39351** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		1/200. Predicted molecular weight: 66 kDa.
IP		Use at an assay dependent dilution.

Target

Function	May play an important role in regulating or promoting cell proliferation in some normal and neoplastically transformed cells.
Pathway	Lipid metabolism; prostaglandin biosynthesis.
Sequence similarities	Belongs to the prostaglandin G/H synthase family. Contains 1 EGF-like domain.
Cellular localization	Microsome membrane. Endoplasmic reticulum membrane.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors