

Product datasheet

Anti-Creatinine antibody ab30719

1 References

Overview

Product name	Anti-Creatinine antibody
Description	Sheep polyclonal to Creatinine
Host species	Sheep
Specificity	Recognizes free creatinine.
Tested applications	Suitable for: Competitive ELISA
Immunogen	Creatinine(N)-BSA

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Storage buffer	pH: 7.20 Preservative: 0.09% Sodium azide Constituents: 0.0536% PBS, 0.87% Sodium chloride
Purity	IgG fraction
Purification notes	Prepared by caprylic acid and ammonium sulphate precipitation techniques.
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab30719** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Competitive ELISA		1/6000.

Target

Relevance	Creatinine is a waste product from protein in the diet and from the muscles of the body. Creatinine is removed from the body by the kidneys; as kidney disease progresses, the level of creatinine in the blood increases.
Cellular localization	Secreted

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors