

Product datasheet

Anti-Cytomegalovirus antibody (HRP) ab87339

Overview

Product name	Anti-Cytomegalovirus antibody (HRP)
Description	Goat polyclonal to Cytomegalovirus (HRP)
Host species	Goat
Conjugation	HRP
Tested applications	Suitable for: WB, ELISA, IHC-P
Species reactivity	Cytomegalovirus AD 169, 100% reactive.
Immunogen	Cytomegalovirus AD 169

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	Preservative: Thimerosal (merthiolate) Constituents: 1% BSA, PBS
Purity	Protein G purified
Purification notes	Purified IgG fraction.
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab87339** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		1/200.
ELISA		1/1000.
IHC-P		Use at an assay dependent dilution.

Target

Relevance

Cytomegalovirus is a member of the herpes virus group, which includes herpes simplex virus types 1 and 2, varicella zoster virus (which causes chicken pox), and Epstein Barr virus (which causes infectious mononucleosis). These viruses share a characteristic ability to remain dormant within the body over a long period.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors