

## Product datasheet

# Donkey F(ab')<sub>2</sub> Anti-Rabbit IgG - H&L, pre-adsorbed ab98486

### Overview

<b>Product name</b>	Donkey F(ab') <sub>2</sub> Anti-Rabbit IgG - H&L, pre-adsorbed
<b>Host species</b>	Donkey
<b>Target species</b>	Rabbit
<b>Tested applications</b>	<b>Suitable for:</b> ICC, ELISA, IHC-P, WB
<b>Minimal cross-reactivity</b>	Chicken, Cow, Goat, Horse, Human, Mouse, Rat, Sheep <a href="#">more details</a>

### Properties

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
<b>Storage buffer</b>	Preservative: 0.09% Sodium Azide Constituents: PBS
<b>Purity</b>	Immunogen affinity purified
<b>Purification notes</b>	This antibody was isolated by affinity chromatography using antigen coupled to agarose beads. F(ab') <sub>2</sub> fragment were generated using a pepsin digestion. Fc fragments and whole IgG molecules have been removed.
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG

### Applications

Our [Abpromise guarantee](#) covers the use of **ab98486** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ICC		Use at an assay dependent dilution.
ELISA		1/100 - 1/30000. Coating: Use 1/100 - 1/500 dilution; Primary: Use 1/1000 - 1/30000 dilution.

Application	Abreviews	Notes
IHC-P		1/200 - 1/2000.
WB		1/1000 - 1/30000. Colorimetric: Use 1/1000 - 1/10000 dilution; Chemiluminescent: Use 1/1000 - 1/30000 dilution.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

### **Terms and conditions**

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors