

Product datasheet

Anti-Duxbl antibody ab105965

1 Image

Overview

| | |
|----------------------------|---|
| Product name | Anti-Duxbl antibody |
| Description | Rabbit polyclonal to Duxbl |
| Host species | Rabbit |
| Tested applications | Suitable for: WB |
| Species reactivity | Reacts with: Mouse |
| Immunogen | Synthetic peptide corresponding to a region within internal sequence amino acids 221-270 (SSQVVQQRSD DQNPNGHLS PTTTPGEQGF HSQPPLQLLT QNRGHNPRES) of Mouse Duxbl (NP_899245). Run BLAST with ExPASy Run BLAST with NCBI |
| Positive control | Mouse lung lysate |

Properties

| | |
|-----------------------------|--|
| Form | Liquid |
| Storage instructions | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. |
| Storage buffer | Preservative: 0.09% Sodium azide Constituents: 2% Sucrose, PBS |
| Purity | Immunogen affinity purified |
| Clonality | Polyclonal |
| Isotype | IgG |

Applications

Our [Abpromise guarantee](#) covers the use of **ab105965** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


| Application | Abreviews | Notes |
|-------------|-----------|--|
| WB | | Use a concentration of 1 µg/ml. Predicted molecular weight: 40 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T. |

Target

Relevance ligand dependent nuclear receptor corepressor-like.

Cellular localization Nuclear

Images



Anti-Duxbl antibody (ab105965) at 1 µg/ml +
Mouse lung lysate at 10 µg

Predicted band size: 40 kDa

Gel concentration: 12%

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors