

Product datasheet

Anti-ENOSF1 antibody ab102474

1 References 1 Image

Overview

Product name	Anti-ENOSF1 antibody
Description	Rabbit polyclonal to ENOSF1
Host species	Rabbit
Tested applications	Suitable for: WB
Species reactivity	Reacts with: Human
Immunogen	Recombinant protein fragment corresponding to a region within amino acids 63-319 of ENOSF1 (AAH01285).
Positive control	A549, HeLa, HepG2 or HCT116 cell lysate

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
Storage buffer	Preservative: 0.01% Thimerosal (merthiolate) Constituents: 20% Glycerol, 0.1M Tris, 0.1M Glycine, pH 7
Purity	Immunogen affinity purified
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab102474** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		1/500 - 1/3000. Predicted molecular weight: 50 kDa.

Target

Function	May regulate thymidylate synthase activity.
Sequence similarities	Belongs to the mandelate racemase/muconate lactonizing enzyme family, ENOSF1 subfamily.
Post-translational modifications	Could be sumoylated.
Cellular localization	Mitochondrion.

Images



Anti-ENOSF1 antibody (ab102474) at 1/1000 dilution + HeLa whole cell lysate at 30 µg

Predicted band size: 50 kDa

10% SDS PAGE

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors