

Product datasheet

Anti-FeLV p27 antibody ab21099

Overview

Product name	Anti-FeLV p27 antibody
Description	Goat polyclonal to FeLV p27
Host species	Goat
Specificity	Non-reactive with normal feline serum.
Tested applications	Suitable for: WB
Immunogen	Purified FeLV p27.
Positive control	Purified FeLV p27 or whole virus to FeLV

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Storage buffer	pH: 7.20 Preservative: 0.1% Sodium azide Constituent: 0.0268% PBS
Purity	Ion Exchange Chromatography
Purification notes	>95% pure. Sodium sulfate precipitation and ion-exchange chromatography.
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab21099** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		Use at an assay dependent dilution. Predicted molecular weight: 60 kDa.

Target

Relevance	The feline leukemia virus (FeLV) is the causative agent of the most important fatal infectious disease complex of domestic cats today. It is an RNA virus belonging to the family Retroviridae.
Cellular localization	Integral membrane protein.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors