

Product datasheet

Goat Anti-Dog IgA H&L ab113729

Overview

Product name	Goat Anti-Dog IgA H&L
Description	Goat polyclonal Secondary Antibody to Dog IgA
Host species	Goat
Target species	Dog
Specificity	By immunoelectrophoresis, the antiserum reacts specifically with Dog IgA. No antibody was detected against immunoglobulin light chains or non-immunoglobulin serum proteins.
Tested applications	Suitable for: Immunoelectrophoresis, Other, Double Immunodiffusion

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
Storage buffer	Preservative: 0.1% Sodium azide Constituent: 99% Whole serum
Purity	Whole antiserum
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab113729** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Immunoelectrophoresis		Use at an assay dependent concentration.
Other		Use at an assay dependent concentration. (ab113729 is suitable for double diffusion).
Double Immunodiffusion		Use at an assay dependent concentration.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors