

Product datasheet

Goat F(ab')₂ Anti-Mouse IgG H&L (SureLight® PE) ab99914

Overview

Product name	Goat F(ab') ₂ Anti-Mouse IgG H&L (SureLight® PE)
Host species	Goat
Target species	Mouse
Tested applications	Suitable for: Flow Cyt
Conjugation	SureLight® PE

Properties

Form	Lyophilised: Reconstitute to 0.5 ml with distilled deionized water, vortex gently and allow vial to sit on ice for 20 minutes.
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	pH: 7.40 Preservative: 0.05% Sodium azide Constituents: 0.1% BSA, 0.87% Sodium chloride, 1.64% Sodium phosphate, 1.71% Sucrose
Purity	IgG fraction
Clonality	Polyclonal
Isotype	IgG
General notes	Lyophilized material is stable for up to one year. After reconstitution should be used within 3 months.

Applications

Our [Abpromise guarantee](#) covers the use of **ab99914** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Flow Cyt		Use a concentration of 3 µg/ml.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors