

Product datasheet

Goat Anti-Mouse IgM mu chain (6nm Gold) preadsorbed ab105292

Overview

Product name	Goat Anti-Mouse IgM mu chain (6nm Gold) preadsorbed
Host species	Goat
Target species	Mouse
Specificity	Recognizes rat IgM. Exhibits minimal crossreactivity to bovine, human and horse serum proteins.
Tested applications	Suitable for: Electron Microscopy, WB, ELISA, Other
Minimal cross-reactivity	Cow, Horse, Human more details
Conjugation	Gold 6nm

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	Preservative: 0.05% Sodium Azide Constituents: 15mg/ml BSA, 0.15M Sodium chloride, 10mM Sodium phosphate, 10mM Sodium borate, pH 8.5
Purity	Immunogen affinity purified
Conjugation notes	6nm Gold
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab105292** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Electron Microscopy		1/10 - 1/20.
WB		Use at an assay dependent dilution.

Application	Abreviews	Notes
ELISA		Use at an assay dependent dilution.
Other		Use at an assay dependent dilution. Light microscopy

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors