


Product datasheet

Anti-Golgin 97 antibody ab112084

[1 Abreviews](#) [1 Image](#)

Overview

Product name	Anti-Golgin 97 antibody
Description	Rabbit polyclonal to Golgin 97
Host species	Rabbit
Tested applications	Suitable for: WB
Species reactivity	Reacts with: Human Predicted to work with: Mouse, Rat, Rabbit, Horse, Chicken, Guinea pig, Cow, Cat, Dog, Zebrafish 
Immunogen	Synthetic peptide, corresponding to N terminal amino acids 18-67 (RPGGATRIPR SVSKESVASM GADSGDDFAS DGSSSREDLS SQLLRNEQI) of Human Golgin 97 (NP_002068). Run BLAST with ExPASy Run BLAST with NCBI
Positive control	HCT15 cell lysate

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Storage buffer	Preservative: 0.09% Sodium azide Constituents: 2% Sucrose, PBS
Purity	>95% by SDS-PAGE
Purification notes	Affinity purifications with biochemical fractionation.
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab112084** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		Use a concentration of 0.2 - 1 µg/ml. Predicted molecular weight: 88 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

Target

Function	Probably involved in maintaining Golgi structure.
Sequence similarities	Contains 1 GRIP domain.
Domain	Extended rod-like protein with coiled-coil domains.
Cellular localization	Golgi apparatus > Golgi stack membrane.

Images



Anti-Golgin 97 antibody (ab112084) at 1 µg/ml
+ HCT15 cell lysate at 10 µg

Predicted band size: 88 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors