

## Product datasheet

# Meis homeobox 3 peptide ab107776

### Overview

---

**Product name** Meis homeobox 3 peptide

### Description

---

**Nature** Synthetic

**Amino Acid Sequence**

**Sequence** PDQNNMWIRDHEDSGSVHLGTPGPSSGGLASQSGDN  
SSDQGDGLDTSVAS

### Specifications

---

Our [Abpromise guarantee](#) covers the use of **ab107776** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

**Applications** Blocking

**Form** Liquid

### Preparation and Storage

---

**Stability and Storage** Shipped at 4°C. Store at -20°C.

### General Info

---

**Sequence similarities** Belongs to the TALE/MEIS homeobox family.  
Contains 1 homeobox DNA-binding domain.

**Cellular localization** Nucleus.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery

- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors