

Product datasheet

Anti-NHLRC2 antibody ab88725

[1 References](#) [2 Images](#)

Overview

Product name	Anti-NHLRC2 antibody
Description	Mouse polyclonal to NHLRC2
Host species	Mouse
Tested applications	Suitable for: WB
Species reactivity	Reacts with: Human
Immunogen	Full length Human protein (AAH32598)
Positive control	Human pancreas lysate and NHLRC2 transfected 293T cell lysate

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Storage buffer	pH: 7.40 Constituent: PBS
Purity	Protein A purified
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab88725** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


Application	Abreviews	Notes
WB		1/500 - 1/1000. Predicted molecular weight: 79 kDa.

Target

Sequence similarities

Contains 6 NHL repeats.
Contains 1 thioredoxin domain.


Images



Anti-NHLRC2 antibody (ab88725) at 1/500 dilution + Human pancreas lysate at 50 µg

Predicted band size: 79 kDa

Observed band size: 45,79 kDa



All lanes : Anti-NHLRC2 antibody (ab88725) at 1/500 dilution

Lane 1 : NHLRC2 transfected 293T cell lysate

Lane 2 : non-transfected 293T cell lysate

Lysates/proteins at 25 µg per lane.

Predicted band size: 79 kDa

Observed band size: 40 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors