

Product datasheet

NMDAR2B peptide ab86148

Overview

**Product name** NMDAR2B peptide

Description

**Nature** Synthetic

**Amino Acid Sequence**

**Sequence** RSQKSPPSIGIAVILV

Specifications

Our [Abpromise guarantee](#) covers the use of **ab86148** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

**Applications** Western blot  
 Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections)  
 Immunocytochemistry  
 Blocking

**Form** Liquid

Preparation and Storage

**Stability and Storage** Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

General Info

**Function** NMDA receptor subtype of glutamate-gated ion channels with high calcium permeability and voltage-dependent sensitivity to magnesium. Mediated by glycine.

**Tissue specificity** Primarily found in the fronto-parieto-temporal cortex and hippocampus pyramidal cells, lower expression in the basal ganglia.

**Sequence similarities** Belongs to the glutamate-gated ion channel (TC 1.A.10.1) family. NR2B/GRIN2B subfamily.

**Cellular localization** Cell membrane. Cell junction > synapse > postsynaptic cell membrane.

## **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors