

Product datasheet

Nitrocellulose Transfer Membrane- 0.22 μ m ab133413

1 References

Overview

Product name	Nitrocellulose Transfer Membrane- 0.22 μ m
General notes	High quality 0.22 μ m nitrocellulose membranes that are more durable than standard nitrocellulose membranes. ab133413 yields clearer blots as the membrane is less prone to distortion. The membranes have been pre-cut for ease of use. Compatible with gels 10 x 10 cm and smaller.

Properties

Form	Solid
Storage instructions	Store at room temperature.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors