

Product datasheet

Rabbit Anti-Cow IgG H&L (Alkaline Phosphatase) ab6928

Overview

Product name	Rabbit Anti-Cow IgG H&L (Alkaline Phosphatase)
Host species	Rabbit
Target species	Cow
Tested applications	Suitable for: Dot blot, ELISA, IHC-P, IHC-Fr, Immunomicroscopy, WB, ICC/IF
Immunogen	Full length native bovine IgG (purified).
Conjugation	Alkaline Phosphatase

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
Storage buffer	pH: 8.00 Preservative: 0.1% Sodium azide Constituents: 1% Polyethylene glycol, 0.00136% Zinc chloride, 0.0095% Magnesium chloride, 0.79% Tris HCl, 50% Glycerol, 0.87% Sodium chloride
Purity	Affinity purified
Purification notes	This product was prepared from monospecific antiserum by immunoaffinity chromatography using bovine IgG coupled to agarose beads
Conjugation notes	Alkaline Phosphatase (Calf Intestine) (Molecular Weight 140,000 daltons)
Clonality	Polyclonal

Applications

Our [Abpromise guarantee](#) covers the use of **ab6928** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Dot blot		

Application	Abreviews	Notes
ELISA		
IHC-P		
IHC-Fr		
Immunomicroscopy		
WB		
ICC/IF		

Application notes

Dot: Use at an assay dependant dilution.

ELISA: 1/1500 - 1/6000.

IHC-P: Use at an assay dependant dilution.

IHC-Fr: Use at an assay dependant dilution.

IM: Use at an assay dependant dilution.

WB: Use at an assay dependant dilution.

Not tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors