

Product datasheet

Rabbit Anti-Rat IgG H&L preadsorbed ab7099

1 References

Overview

Product name	Rabbit Anti-Rat IgG H&L preadsorbed
Host species	Rabbit
Target species	Rat
Tested applications	Suitable for: IHC-Fr, ICC/IF, WB, IHC-P, IP, Immunodiffusion, Conjugation, ELISA
Minimal cross-reactivity	Human more details
Immunogen	Full length native Rat IgG (purified).

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
Storage buffer	Preservative: 0.01% Sodium Azide Constituents: 0.15M Sodium Chloride, 0.02M Potassium Phosphate. pH 7.2
Purity	Affinity purified
Purification notes	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Rat IgG coupled to agarose beads.
Clonality	Polyclonal

Applications

Our [Abpromise guarantee](#) covers the use of **ab7099** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
IHC-Fr		Use at an assay dependent concentration.
ICC/IF		1/15000 - 1/75000.
WB		1/30000 - 1/150000.

Application	Abreviews	Notes
IHC-P		Use at an assay dependent concentration.
IP		Use at an assay dependent concentration.
Immunodiffusion		Use at an assay dependent concentration.
Conjugation		Use at an assay dependent concentration.
ELISA		1/300000 - 1/1.4e+006.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors