

Product datasheet

Rat IgG1, kappa [RTK2071] (FITC) - isotype control ab18404

Overview

Product name	Rat IgG1, kappa [RTK2071] (FITC) - isotype control
Conjugation	FITC. Ex: 493nm, Em: 528nm
Conjugation notes	The antibody was conjugated with FITC under optimal conditions. The solution is free of unconjugated FITC.
Specificity	The RTK2071 immunoglobulin reacts with KLH. The RTK2071 immunoglobulin was chosen as an isotype control after screening on a variety of resting, activated, live, and fixed mouse and human tissues.
Tested applications	Suitable for: Flow Cyt
Immunogen	Trinitrophenol and KLH.

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	Preservative: 0.09% Sodium Azide Constituents: PBS, pH 7.2
Purity	Immunogen affinity purified
Clonality	Monoclonal
Clone number	RTK2071
Isotype	IgG1
Light chain type	kappa

Applications

Our [Abpromise guarantee](#) covers the use of **ab18404** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abbreviations	Notes
Flow Cyt		

Application notes

Flow Cyt: Use at an assay dependent dilution for surface or intracellular staining.

Not tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors