

Product datasheet

Rat monoclonal [LO-MG3-13] Anti-Mouse IgG3 heavy chain ab11656

Overview

Product name	Rat monoclonal [LO-MG3-13] Anti-Mouse IgG3 heavy chain
Host species	Rat
Target species	Mouse
Specificity	ab11656 recognises murine IgG3 (gamma heavy chain) and does not cross react with other murine immunoglobulin classes or subclasses.
Tested applications	Suitable for: ELISA
Immunogen	Polyclonal mouse IgG.

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
Storage buffer	Preservative: 0.09% Sodium Azide Constituents: PBS
Purity	Protein G purified
Purification notes	Purified IgG prepared by affinity chromatography.
Clonality	Monoclonal
Clone number	LO-MG3-13
Myeloma	IR983F
Isotype	IgG1
General notes	This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.

Applications

Our [Abpromise guarantee](#) covers the use of **ab11656** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ELISA		Use a concentration of 5 µg/ml.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors