

Product datasheet

Recombinant Human Bax inhibitor 1 protein ab152735

1 Image

Overview

Product name	Recombinant Human Bax inhibitor 1 protein
Protein length	Full length protein

Description

Nature	Recombinant
Source	Wheat germ

Amino Acid Sequence

Accession	P55061
Species	Human
Sequence	MNIFDRKINF DALLKFSHITPSTQQHLKKVYASFALCMFVAAAGAYVH MV THFIQAGLLSALGSLIMWLMATPHSHETE QKRLGLLAGFAFLTGVGLG PALEFCIAVNPSILPTAFMG TAMI FTCTLSALYARRRSYFLGGILMSA LSLLLLSSLG NVFFGSWLFQANLYVGLVVMCGFVLFDTQLIIEKAEHGD QDYWHCIDLFLDFITVFRKLM MILAMNEKDKKKEKK
Molecular weight	53 kDa including tags
Amino acids	1 to 237

Specifications

Our [Abpromise guarantee](#) covers the use of **ab152735** in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	ELISA SDS-PAGE Western blot
Form	Liquid
Additional notes	Protein concentration is above or equal to 0.05 µg/µl.

Preparation and Storage

Stability and Storage	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.
------------------------------	---


pH: 8.00

Constituents: 0.31% Glutathione, 0.79% Tris HCl

General Info

Function	Suppressor of apoptosis.
Tissue specificity	Highly abundant in testis.
Sequence similarities	Belongs to the B11 family.
Cellular localization	Membrane.

Images



12.5% SDS-PAGE analysis of ab152735 stained with Coomassie Blue.

SDS-PAGE - Bax inhibitor 1 protein (Tagged)
(ab152735)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors