

## Product datasheet

# Recombinant Human CCL27 protein ab93953

1 Image

### Overview

---

**Product name** Recombinant Human CCL27 protein

**Protein length** Full length protein

### Description

---

**Nature** Recombinant

**Source** Escherichia coli

### Amino Acid Sequence

**Species** Human

**Sequence**  
 MFLPPSTAC CTQLYRKPLS DKLLRKVIQV  
 ELQEADGDCH LQAFVLHLAQ RSICIHPQNP  
 SLSQWFEHQE RKLHGTLPKL NFGMLRKMG

### Specifications

---

Our [Abpromise guarantee](#) covers the use of **ab93953** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

**Applications** Mass Spectrometry

SDS-PAGE

**Purity** > 95 % SDS-PAGE.

ab93953 was purified using conventional chromatography techniques.

**Form** Liquid

**Additional notes** Endotoxin Level: < 1.0 EU per 1 µg of protein (determined by LAL method)

### Preparation and Storage

---

**Stability and Storage** Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

Preservative: None

Constituents: 10% Glycerol, 10mM Sodium citrate, pH 3.5

## General Info


---

<b>Function</b>	Chemotactic factor that attracts skin-associated memory T-lymphocytes. May play a role in mediating homing of lymphocytes to cutaneous sites. Binds to CCR10.
<b>Tissue specificity</b>	Testis, thymus, placenta, ovary and skin.
<b>Sequence similarities</b>	Belongs to the intercrine beta (chemokine CC) family.
<b>Cellular localization</b>	Secreted.

---

## Images

---



15% SDS-PAGE analysis of ab93953 (3µg)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors