

Product datasheet

Recombinant Human CITED2 protein ab80152

1 Image

Overview

Product name	Recombinant Human CITED2 protein
Protein length	Protein fragment

Description

Nature	Recombinant
Source	Escherichia coli
Amino Acid Sequence	
Species	Human
Amino acids	1 to 160

Specifications

Our [Abpromise guarantee](#) covers the use of **ab80152** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	Western blot
	ELISA
	SDS-PAGE

Form	Lyophilised
-------------	-------------

Preparation and Storage


Stability and Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. Constituent: PBS
------------------------------	---

Reconstitution	Redissolve the powder with 1ml sterile water will restore to the original concentration 1mg/ml.
-----------------------	---

General Info

Function	Interferes with the binding of transcription factors HIF-1a and STAT2 to p300/CBP.
Sequence similarities	Belongs to the CITED family.

Images



20% SDS PAGE showing ab80152 at approximately 26 kDa.
Coomassie blue staining.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors