

Product datasheet

Recombinant Human CLEC18C protein ab166154

1 Image

Overview

Product name	Recombinant Human CLEC18C protein
Protein length	Full length protein

Description

Nature	Recombinant
Source	Wheat germ

Amino Acid Sequence

Species	Human
Sequence	MAGALNRKESFLLLSLHNRLRSWVQPPAADRRRLDWSDSLAQLAQARAAL CGIPTPSLASGLWRTLQVGWNMQLLPAGLASFVEVSLWFAEGQRYSHAA GECARNATCTHYTQLVWATSSQLGCGRHLC SAGQAAIEAFVCAYSPRGNW EVNGKTMPYKKGAWCSLCTASVSGCFKAWDHAGGLCEVPRNPCRMSCQN HGRLNISTCHCHCPPGYTGRYCQVRC SLQCVHGRFREEECSCVCDIGYGG AQCATKVHFPFHTCDLRIDGDCFMVSSEADTYRARMKCQVTVTPSFLGT P
Amino acids	1 to 301
Tags	proprietary tag N-Terminus

Specifications

Our [Abpromise guarantee](#) covers the use of **ab166154** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	ELISA Western blot
Form	Liquid
Additional notes	Protein concentration is above or equal to 0.05 mg/ml.

Preparation and Storage

Stability and Storage	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.
------------------------------	---

