

Product datasheet

Recombinant Human Eotaxin 3 protein ab87891

1 Image

Overview

---

**Product name** Recombinant Human Eotaxin 3 protein  
**Protein length** Full length protein

Description

---

**Nature** Recombinant  
**Source** Escherichia coli

Amino Acid Sequence

**Species** Human

Specifications

---

Our [Abpromise guarantee](#) covers the use of **ab87891** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

**Applications** SDS-PAGE  
**Purity** > 90 % SDS-PAGE.  
 ab87891 was purified by conventional chromatography, after refolding of the isolated inclusion bodies in a renaturation buffer. Endotoxin level: < 1.0 EU per 1 microgram of protein (determined by LAL method).  
**Form** Liquid

Preparation and Storage

---

**Stability and Storage** Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.  
 Preservative: None  
 Constituents: 10% Glycerol, PBS, pH 7.4

General Info

---


**Function** Chemotactic for eosinophils and basophils. Binds to CCR3.

<b>Tissue specificity</b>	Ubiquitously expressed at low levels in various tissues including heart and ovary.
<b>Sequence similarities</b>	Belongs to the intercrine beta (chemokine CC) family.
<b>Cellular localization</b>	Secreted.

---

## Images

---



15% SDS-PAGE showing ab87891 (3 $\mu$ g).

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors