

Product datasheet

Recombinant Human FAM107B protein ab171700

1 Image

Overview

| | |
|-----------------------|-----------------------------------|
| Product name | Recombinant Human FAM107B protein |
| Protein length | Full length protein |

Description

| | |
|---------------|------------------|
| Nature | Recombinant |
| Source | Escherichia coli |

Amino Acid Sequence

| | |
|-------------------------|---|
| Accession | Q9H098 |
| Species | Human |
| Sequence | MGSSHHHHHH SSGLVPRGSH MGSMAEPDYI EDDNPELIRP QKLINPVKTS RNHQDLHREL LMNQKRGLAP QNKPELQKVM EKRKRDQVIK QKEEEAQKKK SDLEIELLKR QQKLEQLELE KQKLQEEQEN APEFVKVKGN LRRTGQEVAQ AQES |
| Molecular weight | 18 kDa including tags |
| Amino acids | 1 to 131 |
| Tags | His tag N-Terminus |

Specifications

Our [Abpromise guarantee](#) covers the use of **ab171700** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| | |
|--------------------------|--|
| Applications | SDS-PAGE Mass Spectrometry |
| Mass spectrometry | MALDI-TOF |
| Purity | >95% by SDS-PAGE . ab171700 is purified using conventional chromatography techniques. |
| Form | Liquid |

Preparation and Storage

Stability and Storage

Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

pH: 8.00


Constituents: 0.02% DTT, 0.32% Tris HCl, 20% Glycerol, 0.88% Sodium chloride

General Info

Sequence similarities

Belongs to the FAM107 family.

Images



15% SDS-PAGE analysis of ab171700 (3 µg).

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors