

Product datasheet

Recombinant Human KRT34 protein ab114955

1 Image

Overview

Product name	Recombinant Human KRT34 protein
Protein length	Protein fragment

Description

Nature	Recombinant
Source	Wheat germ

Amino Acid Sequence

Accession	O76011
Species	Human
Sequence	<p>SLENTL TESEAHYSSQLSQVQSLITNVESQLAEIRCDLERQNQEYQVLLD</p> <p>VRARLECEINTYRSLLSEEDCKLPCNPCATTNASGNSCGPCGTSQKGCCN</p>
Molecular weight	37 kDa including tags
Amino acids	337 to 436

Specifications

Our [Abpromise guarantee](#) covers the use of **ab114955** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	<p>ELISA</p> <p>SDS-PAGE</p> <p>Western blot</p>
---------------------	--

Form	Liquid
-------------	--------

Additional notes	<p>Protein concentration is above or equal to 0.05 mg/ml.</p> <p>Best used within three months from the date of receipt.</p>
-------------------------	--

Preparation and Storage


Stability and Storage	<p>Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.</p> <p>pH: 8.00</p>
------------------------------	--

Constituents: 0.3% Glutathione, 0.79% Tris HCl

General Info

Tissue specificity	Expressed in the hair follicle.
Sequence similarities	Belongs to the intermediate filament family.

Images



12.5% SDS-PAGE showing ab114955 at approximately 36.63kDa.
Stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors