

Product datasheet

Recombinant Human TMEM175 protein ab164555

1 Image

Overview

Product name Recombinant Human TMEM175 protein
Protein length Full length protein

Description

Nature Recombinant
Source Wheat germ
Amino Acid Sequence
Species Human

Sequence

```

MSQPRTPEQALDTPGDCPPGRRDEDAGEGIQCSQRM
LSFSDALLSIATV
MILPVTHTEISPEQQFDRSVQRLLATRIAVYLMTFLMIV
AWAAHTRLFQ
VVGKTDDTLALLNLACMMTITFLPYTFSLMVTFPDVPL
GIFLFCVCVIAI
GVVQALVGYAFHFPHLLSPQIQRSAHRALYRRHVLGIV
LQGPALCFAAA
IFSLFFVPLSYLLM/TVILLPYVSKVTGWCRDRLLGHRE
PSAHPVEVFSF
DLHEPLSKERVEAFSDGVYAMATLLILDICEDNVPDPK
DVKERFSGSLV
AALSATGPRFLAYFGSFATVGLLWFAHHSFLHVRKAT
RAMGLLNTLSLA
FVGGLPLAYQQTSAFARQPRDELERVRVSTIIFLASIF
QLAMWTTALLH
QAETLQPSVWFGGREHVL MFAKLALYPCASLLAFAS
T
CLLSRFSVGIFHL
MQIAVPCAFLLLRLLVGLALATLRVLRGLARPEHPPPA
PTGQDDPQSPLL PAPC
    
```

Amino acids 1 to 504
Tags proprietary tag N-Terminus

Specifications

Our [Abpromise guarantee](#) covers the use of **ab164555** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	Western blot ELISA
Form	Liquid
Additional notes	Protein concentration is above or equal to 0.05 mg/ml. Previously labelled as MGC4618.


Preparation and Storage

Stability and Storage	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.31% Glutathione, 0.79% Tris HCl
------------------------------	--

General Info

Sequence similarities	Belongs to the TMEM175 family.
Cellular localization	Membrane.

Images



ab164555 on a 12.5% SDS-PAGE stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors