

Product datasheet

Recombinant Human NUP62CL protein (denatured)
ab177611

1 Image

Overview

Product name	Recombinant Human NUP62CL protein (denatured)
Protein length	Full length protein
Description	Recombinant Human NUP62CL protein

Description

Nature	Recombinant
Source	Escherichia coli
Amino Acid Sequence	
Accession	Q9H1M0-2
Species	Human
Sequence	MGSSHHHHHH SSGLVPRGSH MGSMQFTSIS NSLTSTAAIG LSFTTSTTTT ATFTTNTTTT ITSGFTVNQN QLLSRGFENL VPYTSTVRFVYMEK
Molecular weight	10 kDa including tags
Amino acids	1 to 72
Tags	His tag N-Terminus
Additional sequence information	Isoform 2 (AAH17799).

Specifications

Our [Abpromise guarantee](#) covers the use of **ab177611** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE
Purity	> 80 % SDS-PAGE.
Form	Liquid

Preparation and Storage

Stability and Storage

Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

Information available upon request.

General Info


Relevance

NUP62CL belongs to the nucleoporin NSP1/NUP62 family. A nucleoporin is a type of protein which facilitates transport through nuclear pores in the nuclear envelope. There are two named isoforms.

Cellular localization

Cytoplasmic and Nuclear

Images



ab177611 (3ug) on a 15% SDS-PAGE.

SDS-PAGE - Human NUP62CL protein fragment
(ab177611)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors