

Product datasheet

Recombinant Human Pleckstrin protein ab116861

1 Image

Overview

Product name	Recombinant Human Pleckstrin protein
Protein length	Protein fragment

Description

Nature	Recombinant
Source	Wheat germ
Amino Acid Sequence	
Accession	P08567
Species	Human
Sequence	PETIDLGALYLSMKDTEKGIKELNLEKDKKIFNHCFTGN CVIDWLVSNSQ VRNRQEGLMIASSLLNEGYLQPAGDMSKSAVDGTAEN PFLDNPDAFYFP DSGFFCEENS
Molecular weight	38 kDa including tags
Amino acids	121 to 230

Specifications

Our [Abpromise guarantee](#) covers the use of **ab116861** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	ELISA SDS-PAGE Western blot
Form	Liquid
Additional notes	Protein concentration is above or equal to 0.05 mg/ml. Best use within three months from the date of receipt of this protein.

Preparation and Storage

Stability and Storage	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.
------------------------------	---

pH: 8.00

Constituents: 0.3% Glutathione, 0.79% Tris HCl

General Info

Function


Major protein kinase C substrate of platelets.

Sequence similarities

Contains 1 DEP domain.

Contains 2 PH domains.

Images



12.5% SDS-PAGE stained with Coomassie Blue showing ab116861 at approximately 37.73 kDa.

SDS-PAGE - Recombinant Human Pleckstrin protein (ab116861)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors