

Product datasheet

Recombinant Human Spir-1 protein ab163422

1 Image

Overview

Product name	Recombinant Human Spir-1 protein
Protein length	Full length protein

Description

Nature	Recombinant
Source	Wheat germ
Amino Acid Sequence	
Species	Human
Sequence	MRLPSKPYSTLPIFSLGPSALQRGESSMRSEKPSTAH HRPLRSIARFSSK SKSMDKSDEELQFPKELMEDWSTMEVCVDCKKFISEI ISSRRSLVLANK RARLKRKTQSFYMSSPGPSEYCPSERTISEI
Amino acids	1 to 131
Tags	proprietary tag N-Terminus

Specifications

Our [Abpromise guarantee](#) covers the use of **ab163422** in the following tested applications.
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	Western blot ELISA
Form	Liquid
Additional notes	Protein concentration is above or equal to 0.05 mg/ml. Previously labelled as SPIRE1.

Preparation and Storage


Stability and Storage	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00
------------------------------	---

Constituents: 0.31% Glutathione, 0.79% Tris HCl

General Info

Function	Acts as a actin nucleation factor, remains associated with the slow-growing pointed end of the new filament. Involved in vesicle transport processes providing a novel link between actin organization and intracellular transport.
Sequence similarities	Belongs to the spire family. Contains 1 KIND domain. Contains 2 WH2 domains.
Domain	Binds to actin monomers via the WH2 domain. The Spir-box targets binding to intracellular membrane structures.
Cellular localization	Cytoplasm > cytoskeleton. Cytoplasm > perinuclear region. Punctate spots in perinuclear region and cytoplasm, co-localised with Rab11.

Images



ab163422 on a 12.5% SDS-PAGE stained with Coomassie Blue.

SDS-PAGE - Recombinant Human Spir-1 protein
(ab163422)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors