

Product datasheet

Recombinant Human Taf4 protein ab159625

1 Image

Overview

---

**Product name** Recombinant Human Taf4 protein  
**Protein length** Protein fragment

Description

---

**Nature** Recombinant  
**Source** Wheat germ

Amino Acid Sequence

**Species** Human  
**Sequence** LQPPVLSLTQPTQVGVGKQGQPTPLVIQQPPKPGALIRPPQVTLTQTPMV  
 ALRQPHNRIMLTTPQQIQLNPLQVPVVKPAVLPGT  
**Amino acids** 715 to 800  
**Tags** proprietary tag N-Terminus

Specifications

---

Our [Abpromise guarantee](#) covers the use of **ab159625** in the following tested applications.  
 The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

**Applications** ELISA  
 Western blot  
**Form** Liquid  
**Additional notes** Protein concentration is above or equal to 0.05 mg/ml.

Preparation and Storage

---

**Stability and Storage** Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.  
 pH: 8.00  
 Constituents: 0.31% Glutathione, 0.79% Tris HCl

General Info

---


**Relevance**

Taf4 makes part of TFIIID, a multimeric protein complex that plays a central role in mediating promoter responses to various activators and repressors. It potentiates transcriptional activation by the AF-2S of retinoic acid, vitamin D3 and thyroid hormone.

**Cellular localization**

Nuclear

---

**Images**

ab159625 on a 12.5% SDS-PAGE stained with Coomassie Blue.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

---

**Our Abpromise to you: Quality guaranteed and expert technical support**

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

---

**Terms and conditions**

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors