

Product datasheet

Recombinant Human ZNF74 protein ab159851

1 Image

Overview

Product name	Recombinant Human ZNF74 protein
Protein length	Full length protein

Description

Nature	Recombinant
Source	Wheat germ

Amino Acid Sequence

Species	Human
Sequence	<p>MEIPAPEPEKTALSSQDPALSLKENLEDISGWGLPEARSKESVSFKDVAV DFTQEEWGQLDSPQRALYRDVMLENYQNLLALGPPLHKPDVISHLERGEE PWSMQREVPRGPCPEWELKAVPSQQQGICKEEPAQEPIMERPLGGAQAWG RQAGALQRSQAAPWAPAPAMAWDVPVEEFPLRCPLFAQQRVPEGGPLLDT RKNVQATEGRTKAPARLCAGENASTPSEPEKFPQVRRQRGAGAGEGEFVC GECGKAFRQSSSLTLHRRWHSREKAYKCDECGKAFTWSTNLEHRRHTG EKPFFCGECGKAFSCHSSLNVHQRIHTGERPYKCSACEKAFSCSSLLSMH LRVHTGEKPYRCGECGKAFNQRTHLTRHHRIHTGEKPYQCGSCGKAFTCH SSLTVHEKIHSGDKPFKCSDCGKAFNSRSLTLHQRTHTEKPFKCADCG KGFSCHAYLLVHRRHSGEKPFKNECGKAFSSHAYLVHRRHTGEKPF DCSQCWKAFSCHSSLVHQRIHTGEKPYKCSECGRAFSQNHCLIKHQKIH SGEKSFKCEKCGEMFNWSSHLTEHQRLHSEGKPLAIQFNKHLSTYYVPG SLLGAGDAGLRDVPIDALDVAKLLCVPPRAGRNFSLGSKPRN</p>
Amino acids	1 to 644
Tags	proprietary tag N-Terminus

Specifications

Our [Abpromise guarantee](#) covers the use of **ab159851** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	ELISA Western blot
Form	Liquid

Additional notes

Protein concentration is above or equal to 0.05 mg/ml.

Preparation and Storage**Stability and Storage**

Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.

pH: 8.00

Constituents: 0.31% Glutathione, 0.79% Tris HCl

General Info**Function**

May play a role in RNA metabolism.

Tissue specificity

Highly expressed in the fetal brain.

Sequence similarities


Belongs to the krueppel C2H2-type zinc-finger protein family.

Contains 12 C2H2-type zinc fingers.

Contains 1 KRAB domain.

Cellular localization

Nucleus.

Images

ab159851 on a 12.5% SDS-PAGE stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors