

Product datasheet

Recombinant Human ZNF96 protein ab160642

1 Image

Overview

Product name	Recombinant Human ZNF96 protein
Protein length	Full length protein

Description

Nature	Recombinant
Source	Wheat germ
Amino Acid Sequence	
Species	Human

Sequence

MASTWAIQAHMDQDEPLEVKIEEEKYTTRQDWDLRKN
 NTHSREVFRQYFR
 QFCYQETSGPREALSRLRELCHQWLRPETHHTKEQILEL
 LVLEQFLTILPE
 ELQAWVQEQHPESGEEVTVLEDLERELDEPGEQVS
 VHTGEQEMFLQETV
 RLRKEGEPSSMSLQSMKAQPKYESPELESQQEQVLDV
 ETGNEYGNLKQEV
 EEMEPHGKTSSKFENDMSKSARCGETREPEEITEEPS
 ACSREDKQPTCDE
 NGVSLTENS DHTEHQRICPGEESYGDDCGKAFSQH
 SHLIEHQRIHTGDR
 PYKCEECGKA FRGRTVLRHKIIHTGEKPYKNECGKA
 FGRWSALNQHQR
 LHTGEKHYHCNDCGKAFSQKAGLFHHIKIHTRDKPYQC
 TQCNKSFRRSI
 LTQHQQVHTGAKPYECNECGKAFVYNSSLVSHQEIIH
 KEKCYQCKEKGKS
 FSQSGLIQHQRIHTGEKPYKCDVCEKAFIQRSLTEHQ
 RIHTGERPYKCD
 KCGKAFTQRSVLTEHQRIHTGERPYKCDCEGNAFRGIT
 SLIQHQRIHTGE
 KPYQCDECGKA FRQRKKT SYKEILLKNHSEPPQAGVNL
 LLSSLIPEWQSCF RRDL

Amino acids 1 to 604

Tags proprietary tag N-Terminus

Specifications

Our [Abpromise guarantee](#) covers the use of **ab160642** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	Western blot ELISA
Form	Liquid
Additional notes	Protein concentration is above or equal to 0.05 mg/ml.


Preparation and Storage

Stability and Storage	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.31% Glutathione, 0.79% Tris HCl
------------------------------	--

General Info

Relevance	May be involved in transcriptional regulation.
Cellular localization	Nucleus

Images



ab160642 on a 12.5% SDS-PAGE stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors