

Product datasheet

Anti-Streptococcus Group B III antibody ab78838

Overview

Product name	Anti-Streptococcus Group B III antibody
Description	Rabbit polyclonal to Streptococcus Group B III
Host species	Rabbit
Tested applications	Suitable for: Agglutination
Species reactivity	Reacts with Streptococcus Group B III
Immunogen	Tissue/ cell preparation of Streptococcus Group B III (heat inactivated)
Positive control	Antigenic suspension
General notes	Reagents should be allowed to stand at 15°C - 25°C for at least 30 minutes prior to use.

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
Storage buffer	Preservative: 0.08% Sodium azide Constituent: Whole serum
Purity	Whole antiserum
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab78838** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Agglutination		Use at an assay dependent concentration. See attached protocol.

Target

Relevance

Streptococcus is a genus of spherical, Gram-positive bacteria of the phylum Firmicutes. Streptococcus agalactiae is a gram-positive streptococcus characterized by the presence of Group B Lancefield antigen, and so takes the name Group B Streptococcus. Infection with Group B Streptococcus (GBS), also known as Streptococcus agalactiae and more colloquially as Strep B and group B Strep, can cause serious illness and sometimes death, especially in newborn infants, the elderly, and patients with compromised immune systems. Group B streptococci are also prominent veterinary pathogens, because they can cause bovine mastitis (inflammation of the udder) in dairy cows. The species name "agalactiae" meaning "no milk", alludes to this.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors