

Product datasheet

Synaptogyrin 3 peptide (131-144) ab45657

Overview

---

**Product name** Synaptogyrin 3 peptide (131-144)

Description

---

**Nature** Synthetic

**Amino Acid Sequence**

**Species** Human

**Sequence** CQRTAPGPATTQAGD

**Amino acids** 131 to 144

Specifications

---

Our [Abpromise guarantee](#) covers the use of **ab45657** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

**Applications** Blocking - Blocking peptide for Anti-Synaptogyrin 3 antibody ([ab26059](#))

**Form** Liquid

Preparation and Storage

---

**Stability and Storage** Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

General Info

---

**Function** Involved in the positive regulation of dopamine transporter activity. Probably facilitates the physical and functional interactions of the transporter with the dopamine vesicular storage system, allowing a more efficient loading of the vesicles with extracellular dopamine after release.

**Tissue specificity** Expressed in brain and placenta.

**Sequence similarities** Belongs to the synaptogyrin family.  
Contains 1 MARVEL domain.

**Cellular localization** Membrane. Cytoplasmic vesicle > secretory vesicle > synaptic vesicle. Cell junction > synapse.  
Found at the neuromuscular synapses.

---

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors