

Product datasheet

Anti-TINP1 antibody ab106484

[3 Images](#)

Overview

Product name	Anti-TINP1 antibody
Description	Chicken polyclonal to TINP1
Host species	Chicken
Specificity	ab106484 is predicted not to cross-react with NSA1.
Tested applications	Suitable for: IHC-P, WB
Species reactivity	Reacts with: Mouse, Rat, Human
Immunogen	A 13 amino acid peptide near the carboxy terminus of Human TINP1 (NP_055701).
Positive control	Rat brain tissue lysate Human brain tissue

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at 4°C (stable for up to 12 months).
Storage buffer	Preservative: 0.02% Sodium Azide Constituents: PBS
Purity	Immunogen affinity purified
Clonality	Polyclonal
Isotype	IgY

Applications

Our [Abpromise guarantee](#) covers the use of **ab106484** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
IHC-P		Use a concentration of 5 µg/ml.
WB		Use a concentration of 1 - 2 µg/ml. Predicted molecular weight: 30 kDa.

Target


Function	Involved in the biogenesis of the 60S ribosomal subunit. May play a part in the quality control of pre-60S particles.
Sequence similarities	Belongs to the ribosomal protein S8e family.
Cellular localization	Nucleus > nucleolus.

Images



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-TINP1 antibody (ab106484)


Immunohistochemistry of TINP1 in paraffin embedded Human brain tissue with ab106484 at 5 µg/ml.



Western blot - Anti-TINP1 antibody (ab106484)

Anti-TINP1 antibody (ab106484) at 2 µg/ml +
Rat brain tissue lysate at 15 µg

Predicted band size: 30 kDa



Western blot - Anti-TINP1 antibody (ab106484)

Lane 1 : Anti-TINP1 antibody (ab106484) at 1 µg/ml

Lane 2 : Anti-TINP1 antibody (ab106484) at 2 µg/ml

All lanes : Rat Brain Tissue Lysate

Lysates/proteins at 15 µg per lane.

Predicted band size: 30 kDa

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors