

Product datasheet

Anti-TIS11B antibody ab35801

[1 References](#) [1 Image](#)

Overview

Product name	Anti-TIS11B antibody
Description	Rabbit polyclonal to TIS11B
Host species	Rabbit
Tested applications	Suitable for: WB, ELISA
Species reactivity	Reacts with: Human Predicted to work with: Mouse, Rat, Cow, Dog
Immunogen	Synthetic peptide corresponding to Human TIS11B aa 66-115. Sequence: QLLSSLKGEPAPALSSRDSRFRDRSFSEGGERLLPTQ KQPGGGQVNSSRY Database link: Q07352 Run BLAST with Run BLAST with
Positive control	Human placenta lysate.

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Storage buffer	Preservative: 0.09% Sodium azide Constituents: 2% Sucrose, PBS
Purity	Immunogen affinity purified
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab35801** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		Use a concentration of 0.5 µg/ml. Predicted molecular weight: 38 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.
ELISA		1/62500.

Target

Function	Probable regulatory protein involved in regulating the response to growth factors.
Sequence similarities	Contains 2 C3H1-type zinc fingers.
Cellular localization	Cytoplasm. Nucleus. Shuttles between the nucleus and the cytoplasm.

Images



Anti-TIS11B antibody (ab35801) at 0.5 µg/ml (in 5% skim milk / PBS buffer) + Human placental lysate at 10 µg

Secondary

HRP-conjugated anti-Rabbit IgG

Predicted band size: 38 kDa

Observed band size: 45 kDa

[why is the actual band size different from the predicted?](#)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors