


Product datasheet

Anti-ZCCHC7 antibody ab104503

[2 Images](#)

Overview

Product name	Anti-ZCCHC7 antibody
Description	Rabbit polyclonal to ZCCHC7
Host species	Rabbit
Tested applications	Suitable for: WB, IHC-P
Species reactivity	Reacts with: Human Predicted to work with: Mouse, Rat 
Immunogen	Recombinant fragment, corresponding to amino acids 83-332 of Human ZCCHC7 (BC036940).
Positive control	293T cell lysate; Human Fetal Brain lysate; Human small intestine

Properties

Form	Lyophilised:Reconstitute in 200ul sterile H2O.
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
Storage buffer	Preservative: 0.02% Sodium Azide Constituents: PBS, pH 7.2
Purity	Protein A purified
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab104503** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


Application	Abreviews	Notes
WB		1/500 - 1/1000. Predicted molecular weight: 63 kDa.
IHC-P		1/100 - 1/500.

Target

Sequence similarities Contains 4 CCHC-type zinc fingers.

Cellular localization Nucleus, nucleolus.

Images




Western blot - Anti-ZCCHC7 antibody (ab104503)

All lanes : Anti-ZCCHC7 antibody (ab104503) at 1/500 dilution

Lane 1 : 293T cell lysate

Lane 2 : Human fetal brain lysate

Predicted band size: 63 kDa



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-ZCCHC7 antibody (ab104503)

ab104503, at 1/100 dilution, staining ZCCHC7 in formalin-fixed, paraffin-embedded Human small intestine by Immunohistochemistry.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors