

Product datasheet

pQM-CTag/B vector ab1058

1 Image

Overview

Product name pQM-CTag/B vector

General notes The CMV promoter is covered under U. S. Patents 5,168,062 and 5,385,839 and its use is permitted for research purpose only. Any other use of the CMV promoter requires a licence from the University of Iowa Research Foundation, 214 Technology Innovation Center, Iowa City, IA 52242.


For more information on the E2 tag system see the datasheets for [ab978](#), [ab977](#), [ab997](#) and [ab1013](#).

Properties

Storage instructions Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Vector notes For more information on the E2 tag system see the datasheets for [ab978](#), [ab977](#), [ab997](#) and [ab1013](#).

Images



pQM-CTag/B map

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors