

Product datasheet

Sheep Anti-Dog IgG H&L (FITC) ab112839

Overview

Product name	Sheep Anti-Dog IgG H&L (FITC)
Host species	Sheep
Target species	Dog
Specificity	No antibody was detected against immunoglobulin light chains or non-immunoglobulin serum proteins.
Tested applications	Suitable for: ICC/IF, IHC-P
Conjugation	FITC. Ex: 493nm, Em: 528nm

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	Preservative: 0.09% Sodium azide Constituents: 99% PBS, 0.2% BSA
Purity	Immunogen affinity purified
Purification notes	Antiserum was solid phase adsorbed to ensure class specificity. ab112839 was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to fluorescein isothiocyanate (FITC).
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab112839** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ICC/IF		1/200 - 1/2000.
IHC-P		1/200 - 1/2000.

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors