

Product datasheet

Anti-YIF1B antibody ab112951

1 Image

Overview

| | |
|----------------------------|---|
| Product name | Anti-YIF1B antibody |
| Description | Rabbit polyclonal to YIF1B |
| Host species | Rabbit |
| Tested applications | Suitable for: WB |
| Species reactivity | Reacts with: Human |
| Immunogen | Synthetic peptide conjugated to KLH, corresponding to a region within N terminal amino acids 81-110 of Human YIF1B (NP_001138935.1, NP_001138934.1, NP_001138933.1, NP_001034762.1, NP_291035.1, NP_001034761.1, NP_001034760.1). |
| Positive control | MCF7 cell line lysate. |

Properties

| | |
|-----------------------------|---|
| Form | Liquid |
| Storage instructions | Shipped at 4°C. Store at 4°C (up to 6 months). Store at -20°C long term. |
| Storage buffer | Preservative: 0.09% Sodium azide Constituent: 99% PBS |
| Purity | Immunogen affinity purified |
| Purification notes | ab112951 is purified through a protein A column, followed by peptide affinity purification. |
| Clonality | Polyclonal |
| Isotype | IgG |

Applications

Our [Abpromise guarantee](#) covers the use of **ab112951** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


| Application | Abreviews | Notes |
|-------------|-----------|--|
| WB | | 1/100 - 1/500. Predicted molecular weight: 34 kDa. |

Target

Sequence similarities Belongs to the YIF1 family.

Cellular localization Membrane.

Images



Anti-YIF1B antibody (ab112951) at 1/100 dilution + MCF7 cell line lysate at 35 µg

Predicted band size: 34 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors