

Product datasheet

ab113724

Overview

Product name

Description

Sheep polyclonal Secondary Antibody to Cow Whole Serum

Host species

Sheep

Target species

Cow

Specificity

ab113724 reacts with multiple serum proteins in Cow serum.

Tested applications

Suitable for: Immunoelectrophoresis, Other, Double Immunodiffusion

Properties

Form

Liquid

Storage instructions

Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

Storage buffer

Preservative: 0.1% Sodium azide
Constituent: 99% Whole serum

Purity

Whole antiserum

Clonality

Polyclonal

Isotype

IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab113724** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Immunoelectrophoresis		Use at an assay dependent concentration.
Other		Use at an assay dependent concentration. (ab113724 is suitable for double diffusion Oucterlony).
Double Immunodiffusion		Use at an assay dependent concentration.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors