

Product datasheet

Anti-FeSODB1 antibody ab115977

Overview

| | |
|----------------------------|--|
| Product name | Anti-FeSODB1 antibody |
| Description | Chicken polyclonal to FeSODB1 |
| Host species | Chicken |
| Tested applications | Suitable for: WB, ELISA |
| Species reactivity | Reacts with: Leishmania |
| Immunogen | Synthetic peptide (Leishmania) Purified from FeSODB1 protein, covalently modified to carrier protein and emulsified in adjuvant. |

Properties

| | |
|-----------------------------|---|
| Form | Liquid |
| Storage instructions | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. |
| Storage buffer | Preservative: 0.02% Sodium azide |
| Purity | Immunogen affinity purified |
| Clonality | Polyclonal |
| Isotype | IgY |

Applications

Our [Abpromise guarantee](#) covers the use of **ab115977** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| Application | Abreviews | Notes |
|-------------|-----------|--|
| WB | | 1/500. Predicted molecular weight: 22 kDa. |
| ELISA | | 1/500. |

Target

| | |
|------------------|--|
| Relevance | Destroys radicals which are normally produced within the cells and which are toxic to biological |
|------------------|--|

systems.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors