

Product datasheet

Anti-NKAIN2 antibody ab126547

1 Image

Overview

<b>Product name</b>	Anti-NKAIN2 antibody
<b>Description</b>	Rabbit polyclonal to NKAIN2
<b>Host species</b>	Rabbit
<b>Tested applications</b>	<b>Suitable for:</b> IHC-P
<b>Species reactivity</b>	<b>Reacts with:</b> Human
<b>Immunogen</b>	antigen sequence: GCTVTSVTPA PDWAPEDHRY ITVSGCLLEY QYIEVAH, corresponding to amino acids 114-150 of Human NKAIN2. <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>
<b>Positive control</b>	Human placenta tissue.

Properties

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
<b>Storage buffer</b>	pH: 7,20 Preservative: 0.02% Sodium azide Constituents: 59% PBS, 40% Glycerol
<b>Purity</b>	Immunogen affinity purified
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab126547** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
IHC-P		1/20 - 1/50. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

## Target


---

<b>Tissue specificity</b>	Expressed in fetal brain. Weakly expressed in adult brain and thymus. Not expressed in any other normal tissue examined.
<b>Involvement in disease</b>	Note=A chromosomal aberration involving NKAIN2 is a cause of lymphoma. Deletion del(6)(q13q21) within NKAIN2 and involving SUSP1 generates the SUSP1-NKAIN2 product.
<b>Sequence similarities</b>	Belongs to the NKAIN family.
<b>Cellular localization</b>	Cell membrane.

---

## Images

---



ab126547 at 1/20 dilution staining NKAIN2 in paraffin embedded Human placenta tissue by Immunohistochemistry.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-NKAIN2 antibody (ab126547)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors