

Product datasheet

Anti-L3HYPDH antibody - C-terminal ab177320

1 Image

Overview

Product name	Anti-L3HYPDH antibody - C-terminal
Description	Rabbit polyclonal to L3HYPDH - C-terminal
Host species	Rabbit
Tested applications	Suitable for: WB
Species reactivity	Reacts with: Human Predicted to work with: Mouse, Rat, Rabbit, Horse, Guinea pig, Cow, Cat, Dog
Immunogen	Synthetic peptide within Human L3HYPDH aa 297-346 (C terminal). The exact sequence is proprietary. (NP_653182). Sequence: AFKSSATGSVFTGKAVREAKCGDFKAVIMEVSGQAHYGTASFIEDDDP Database link: Q96EM0 Run BLAST with Run BLAST with
Positive control	OVCAR-3 whole cell lysate.

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
Storage buffer	Constituents: 2% Sucrose, 98% PBS
Purity	Immunogen affinity purified
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab177320** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
-------------	-----------	-------

WB Use a concentration of 1 µg/ml. Predicted molecular weight: 38 kDa.

Target

Images



Anti-L3HYPDH antibody - C-terminal (ab177320) at 1 µg/ml + OVCAR-3 whole cell lysate at 10 µg

Predicted band size: 38 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors