

Product datasheet

Anti-HTLV I p19 antibody [6E10/B8] ab18138

Overview

| | |
|----------------------------|--|
| Product name | Anti-HTLV I p19 antibody [6E10/B8] |
| Description | Mouse monoclonal [6E10/B8] to HTLV I p19 |
| Host species | Mouse |
| Specificity | The clone number has been updated from (8.A.024) to (6E10/B8) both clone numbers name the same antibody clone. |
| Tested applications | Suitable for: ELISA, WB |
| Immunogen | HTLV1 gag p19, Recombinant (Human T-Cell Leukemia Virus). |

Properties

| | |
|-----------------------------|--|
| Form | Liquid |
| Storage instructions | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Add glycerol to a final volume of 50% for extra stability and aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle. |
| Storage buffer | Constituents: PBS, 0.038% EGTA, 0.0292% EDTA |
| Purity | IgG fraction |
| Purification notes | = 95% Purified by ion-exchange chromatography and gel filtration. |
| Clonality | Monoclonal |
| Clone number | 6E10/B8 |
| Isotype | IgG |

Applications

Our [Abpromise guarantee](#) covers the use of **ab18138** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| Application | Abreviews | Notes |
|-------------|-----------|-------------------------------------|
| ELISA | | Use at an assay dependent dilution. |
| WB | | 1/50 - 1/500. |

Target

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors